BRAEBON MEDICAL CORPORATION	
Unites States of America: 102-2981 Ford Street Extension Ogdensburg, NY 13669-3474 1.888.462.4841	Canada: 1-100 Schneider Road, Kanata, ON, Canada, K2K 1Y2 01.613.831.6690
AUTHORIZED EUROPEAN REPRESENTATIVE	
CE	CEpartner4U BV Estdooriaan 13 3951 DB Mearn The Netherlands Tel: +31.343.442.524
© 2008 BRAEBON® Medical Corporation	

en	English
Snore Sensor and Snore Microphone User Guide	
All Models	

Intended Use

The BRAEBON Medical Corporation *Ultima* Snore Sensor™and Snore Microphone are intended for use during sleep disorder studies as a measure of snoring for recording onto a data acquisition system. The snore sensor (red disc) uses piezo-electric technology and plugs directly into the electrode headbox; the snore microphone (red hourglass shape) uses a condenser electret microphone and plugs directly into the electrode headbox.

Warning: Snore sensors and microphones are for diagnostic purposes only and are not intended as apnea (apnoea) monitors.

Caution: US Federal law restricts this device to sale by or on the order of a physician. Never steam autoclave the devices. Do not immerse the snore microphone in any liquid. Use isopropyl alcohol pads to gently clean the snore sensor and snore microphone.

Note: Refer to www.braebon.com for more information. Product specifications are subject to change without notice.

How to Use

Refer to your sleep system user guide for manufacturer’s recommended settings. Red 1.5 mm touchproof connector is positive; black is negative.

General recording system settings:

- Low frequency filter: 10 Hz
- High frequency filter: 50 Hz to 100 Hz
- 50 / 60 Hz notch filter: optional

- Sensitivity / Gain: 27.5 uV/mm or 15,000 gain (increase or decrease as needed) Attach the snore microphone on the nasion, cheek or chin using medical grade tape. Ensure the snore microphone is pointing towards the nose and mouth. Attach the snore sensor (circular shape) to the side of the neck at point of maximal snoring vibration. Ensure the nipple side is face down against skin.

Specifications

Operating Temperature: 5°C (41°F)-32°C (90°F)

Storage Temperature: 0°C (32°F)-60°C (140°F)

Operating/Storage Humidity: 10–95% Non-condensing

Disinfection and Cleaning

Caution: Never steam autoclave the devices. Do not immerse the snore microphone in water. Use isopropyl alcohol pads to clean the sensors.

The *Ultima* Snore Sensor™ and *Ultima* Snore Microphone™ are classified as non-critical items because they only contact intact skin. You may use isopropyl alcohol pads to gently clean the snore sensor and snore microphone. Do not immerse the snore microphone in any liquid. For additional information, refer to the APIC guidelines for selection and use of disinfectants (American Journal of Infection Control, Vol. 18, number 2, April, 1990, p. 99-117).

Warranty

BRAEBON MEDICAL CORPORATION warrants to the first consumer that the snore sensor (the “Sensor”), when shipped in its original container, will be free from defective workmanship, performance and materials and agrees that it will, at its option, either repair the defect or replace the defective Sensor or part thereof at no charge to the purchaser for parts or labor for a time period of nine months from the date of purchase; the snore microphone (the “Sensor”) has a warranty of six months from date of purchase. The warranty described herein shall be the sole and exclusive warranty granted by BRAEBON MEDICAL CORPORATION and shall be the sole and exclusive remedy available to the purchaser. Use of the Sensor constitutes total and complete acceptance of this warranty. Correction of defects, in the manner and for the time period described herein, shall constitute complete fulfillment of all liabilities and responsibilities of BRAEBON MEDICAL CORPORATION to the purchaser with respect to the Sensor and shall constitute full satisfaction of all claims, whether based on contract, negligence, strict liability or otherwise. In no event shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any loss of revenues or damage, direct, incidental, or consequential, including property damage, loss of profit, or personal injury resulting from the use or misuse of, or the inability to use this product. Nor shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any damages or defects in the Sensor which were caused by abuse, misuse, tampering, neglect, incorrect battery type, or repairs or attempted repairs performed by anyone other than an authorized servicer.

de	Deutsch

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