<p></p>
<p>en English</p>
<p>Q-RIP™ Respiratory Effort System User Guide</p>
<p>All Models</p>

Intended Use

The BRAEBON MEDICAL CORPORATION Q-RIP™ Respiratory Effort System™ is intended for use during sleep disorder studies as a measure of respiratory effort for recording onto a data acquisition system. The sensor uses respiratory inductance plethysmography technology and plugs directly into the electrode headbox or a DC amplifier.

Warning: Respiratory effort sensors are for diagnostic purposes only and are not intended as an apnea (apnoea) monitor.

Caution: US Federal law restricts this device to sale by or on the order of a physician. You may wash the belt, but not the interface.

Note: Refer to www.braebon.com for more information. Product specifications are subject to change without notice.

How to Use

Refer to your sleep system user guide for manufacturer's recommended settings. Red 1.5 mm touchproof connector is positive; black is negative. General recording system settings:

- Low frequency filter: 0.1 Hz or lower (TC = 1 sec or longer)
- High frequency filter: 5 Hz to 15 Hz
- 50 / 60 Hz notch filter: optional
- Sensitivity / Gain: 2 mV/cm or 4,000 gain (increase or decrease as needed)

1) Q-RIP Belt Placement: Adjust the Q-RIP™ Respiratory Effort Sensor™ to fit the patient comfortably. Place the belt snugly around either the patient's chest or abdomen. The chest belt should be just below or above the nipple line and the abdominal belt should be below the rib cage. To use the buckle connector, simply fasten one end of the belt into the other end of the belt. Ensure that the belt is neither excessively tight nor loose.

2) Using the Q-RIP™ Interface: To use the Q-RIP™ interface, turn the unit on and plug the two safety pin connectors into the chest and abdomen belts. The connectors may be found next to the buckles. Next, plug the abdomen belt into the interface abdomen input and the chest belt into the interace chest input. Connect the abdomen and chest outputs to your sleep recording system. Use low gain if plugging into a headbox and amplifier; use high gain for a DC connection. You may or may not choose to use the sum output.

Two AA alkaline batteries will last approximately 70 8-hour recordings (i.e., 70 nights). The interface will automatically turn off after ten hours of continuous use. The LED will blink red when it is time to change the batteries.

Adjust the belts and amplifier sensitivity / gain as required.

Specifications

Operating Temperature: 5°C (41°F)-32°C (90°F)
Storage Temperature: 0°C (32°F)-60°C (140°F)
Operating/Storage Humidity: 10–95% Non-condensing

Disinfection and cleaning

Caution: US Federal law restricts this device to sale by or on the order of a physician. You may wash the belt, but not the interface.

The *Ultima* Q-RIP Respiratory Effort System™ is classified as a non-critical item because it only contacts clothing and intact skin. You may wash the respiratory belt, but not the interface. For additional information, refer to the APIC guidelines for selection and use of disinfectants (American Journal of Infection Control, Vol. 18, number 2, April, 1990, p. 99-117).

Ensure that The Q-RIP Respiratory Effort Sensor™ is completely dry before using. Never immerse the Q-RIP™ interface into any liquids at any time.

Steam autoclaving is not recommended.

Warranty

BRAEBON MEDICAL CORPORATION warrants to the first consumer that the Q-RIP™ respiratory effort system (the "Sensor"), when shipped in its original container, will be free from defective workmanship, performance and materials and agrees that it will, at its option, either repair the defect or replace the defective Sensor or part thereof at no charge to the purchaser for parts or labor for a time period of one year from the date of purchase. The warranty described herein shall be the sole and exclusive warranty granted by BRAEBON MEDICAL CORPORATION and shall be the sole and exclusive remedy available to the purchaser. Use of the Sensor constitutes total and complete acceptance of this warranty. Correction of defects, in the manner and for the time period described herein, shall constitute complete fulfillment of all liabilities and responsibilities of BRAEBON MEDICAL CORPORATION to the purchaser with respect to the Sensor and shall constitute full satisfaction of all claims, whether based on contract, negligence, strict liability or otherwise. In no event shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any loss of revenues or damage, direct, incidental, or consequential, including property damage, loss of profit, or personal injury resulting from the use or misuse of, or the inability to use this product. Nor shall BRAEBON MEDICAL CORPORATION be liable, or in any way responsible, for any damages or defects in the Sensor which were caused by abuse, misuse, tampering, neglect, incorrect battery type, or repairs or attempted repairs performed by anyone other than an authorized servicer.

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